


## TRANSFER ANDENTREPRENEUR CENTRE

## **Confocal microscope**

The confocal microscope allowes users to analyze the surface structure of various components.

Training period for independent usage: approx. 8 h

Nanofocus µsurf expert



The Nanofokus µsurf expert is a confocal microscope with which surface structures can be visualized.

Magnification
Camera resolution

Software

Maximum Z-Hub

10x; 20x; 50x

4 megapixel at 16 bit color depth (HDR-function) Installed on control computer: μsoft metrology, μsoft analysis

150 mm

TRANSFER ANDENTREPRENEUR CENTRE